


## Radiation Monitoring System DLMon

### Fields of Application:

- Dose rate monitoring system for use in nuclear medicine
- Area monitoring ( e.g. manufacturing of radionuclides, laboratories, nuclear power plants, or storage facilities of radioactive materials)

### Features:


- Custom-oriented measuring system consisting of intelligent radiation detectors and displays (up to 255 single elements).
- Data transfer with RS485 bus system
- Analysis software based on advanced Client-Server-Concept (custom-designed analysis of acquired data)
- Data storage with Microsoft SQL-SERVER
- Customization of software



**Address:** Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld


**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)



<p><b><u>Intelligent radiation Detector:</u></b></p> <p><b><u>SON16</u></b></p>	
Measuring value:	Photon dose equivalent rate or ambient dose equivalent rate
Detector:	Geiger-Muller tube Typ 70004
Measurement range:	1 $\mu\text{Sv/h}^*$ ... 500 $\mu\text{Sv/h}$ * (Measurement time 200 s, Measurement uncertainty < 3 %)
Energy range:	35 keV ... 1.3 MeV
Collimator:	Material: Lead Acceptance angle: $\pm 30^\circ$
Electronics:	<ul style="list-style-type: none"> <li>• micro-controller with internal memory</li> <li>• 128-kByte-Flash-EEPROM</li> <li>• data transfer by RS-485 interface</li> <li>• calculation of dose rate from counting rate</li> <li>• correction for detector dead time and background counting rate</li> <li>• data memory (up to 9088) values including date and time in a non-volatile cyclic memory</li> <li>• acquisition parameters (acquisition time, number of cycles) can be defined</li> <li>• continuous self-monitoring of all internal and external operating voltages;</li> <li>• functional tests of the detectors, interfaces and memories</li> </ul>

**Address:** STEP Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld


**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)

<p><b><u>Intelligent radiation Detector:</u></b></p> <p><b><u>SON20</u></b></p>	
Measuring value:	Photon dose equivalent rate or ambient dose equivalent rate
Detectors:	Geiger-Muller tube Typ 70013E (*)
Measurement range:	0.5 $\mu\text{Sv/h}^*$ ... 200 $\mu\text{Sv/h}$ * (Measurement time 200 s, Measurement uncertainty < 3 %)
Energy range:	35 KeV ... 1.3 MeV
Collimator:	none
Electronics:	<ul style="list-style-type: none"> <li>• micro-controller with internal memory</li> <li>• 128-kByte-Flash-EPROM</li> <li>• data transfer by RS-485 interface</li> <li>• calculation of dose rate from counting rate</li> <li>• correction for detector dead time and background counting rate</li> <li>• data memory (up to 9088) values including date and time in a non-volatile cyclic memory</li> <li>• acquisition parameters (acquisition time, number of cycles) can be defined</li> <li>• continuous self-monitoring of all internal and external operating voltages;</li> <li>• functional tests of the detectors, interfaces and memories</li> </ul>

(\*) The Type SON20 is available in versions with different GM-Tubes and several measurement ranges.



**Address:** STEP Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld

**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)

<p><b><u>Intelligent radiation Detector:</u></b></p> <p><b><u>SON31</u></b></p>	
Measuring value:	Ambient dose equivalent rate H*(10)
Detectors:	Geiger-Muller tube Typ 70031A
Measurement range:	0.1 $\mu\text{Sv/h}^*$ ... 1000 $\mu\text{Sv/h}$ * (Measurement time 200 s, Measurement uncertainty < 3 %)
Energy range:	35 KeV ... 1.3 MeV
Collimator:	optional
Degree of protection:	IP 65 (die-cast case / weatherproof protection)
Electronics:	<ul style="list-style-type: none"> <li>• micro-controller with internal memory</li> <li>• 128-kByte-Flash-EEPROM</li> <li>• data transfer by RS-485 interface</li> <li>• calculation of dose rate from counting rate</li> <li>• correction for detector dead time and background counting rate</li> <li>• data memory (up to 9088) values including date and time in a non-volatile cyclic memory</li> <li>• acquisition parameters (acquisition time, number of cycles) can be defined</li> <li>• continuous self-monitoring of all internal and external operating voltages;</li> <li>• functional tests of the detectors, interfaces and memories</li> </ul>

**Address:** STEP Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld

**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)

<p><b><u>Intelligent display unit:</u></b></p> <p>AT1 – 1</p> <p>AT1 - 2</p>	
<p><b><u>Intelligent big display unit:</u></b></p> <p>GAT-1</p>	
<p>Single display for each detector or double display for two detectors, resp.</p>	
<p>Display of the measuring value:</p> <ul style="list-style-type: none"> <li>• 3 1/2 - digit display (height of digits: 17.8 mm), LCD with backlighting</li> <li>• quasi-analogue bar graph (20 elements)</li> <li>• optical and acoustic error signalling (reset by software)</li> <li>• 3 LED (red, yellow, green) for rapid signalling of the health risk for staff (AT1-x)</li> <li>• 3 different colors for rapid signalling of the health risk for staff (GAT1)</li> <li>• data transfer by RS-485 interface</li> <li>• EPROM for data storage</li> </ul>	
<p>Electronics:</p>	<ul style="list-style-type: none"> <li>• micro-controller with internal memory</li> <li>• 128-kByte-Flash-EPROM</li> <li>• data transfer by RS-485 interface</li> <li>• data memory (up to 9088) values including date and time in a non-volatile cyclic memory</li> <li>• internal functional tests</li> </ul>



**Address:** STEP Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld

**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)

## **RMon-client -software:**

Analysis software for room monitoring for use in laboratories and exposed rooms.

- Windows based software for PC (Window 2000, XP, 7)
- Display of the arrangement of each detector in the rooms / process to be monitored (customized ground view / floor plan)
- Detector overview: display current dose rate of all detectors as bar chart as well as by a colour code
- Display the diagram with dose rate vs. time for each detector
- Calculate the dose
- Detector calibration
- Marginal check: if pre-selected thresholds are exceeded-emission of warning signals and/or activation of safety facilities





**Address:** STEP Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld

**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)

## **DLMon-client-software:**

Analysis software for patient monitoring for use in radioiodine therapy.

- Display of the arrangement of each detector in the rooms to be monitored
- Detector overview: display of the current dose rate digitally as well as by a colour code
- Dose rate histogram
- Long-term archiving of acquired data
- Detector calibration
- Automatic monitoring of all components and functions
- Marginal check: if pre-selected thresholds are exceeded - emission of warning signals and/or activation of safety facilities
- Patient check-in (optional: chip-card reader)
- Patient data base
- Diagrams including the display of discharge-thresholds (user programmable)
- Data fit to estimate the effective half-life of the incorporated nuclide
- Transfer of patients by mouse-click
- Printing forms for patient discharge
- Drawing up reports on chosen periods



**Address:** STEP Sensortechnik und  
Elektronik Pockau GmbH  
Siedlungsstrasse 5-7  
D-09509 Pockau-Lengefeld

**Phone:** (0)37367 / 9791 and 9792  
**Fax:** (0)37367 / 77730  
**E-Mail:** [info@step-sensor.de](mailto:info@step-sensor.de)  
**URL:** [www.step-sensor.de](http://www.step-sensor.de)